



## DRY BATH INCUBATOR BDIB-202

# DRY BATH INCUBATOR BDIB-202

## DRY BATH INCUBATOR

Dry bath incubators are ideal for incubation and activation of cultures, enzyme reactions, immunoassays, melting / boiling points and a wide variety of other laboratory procedures.



Precise temperature control, fast and safe heating, temperature deviation calibration function.

Multiple block options, support block customization.

Convenient block replacement, easy to clean and sterilize.

## SPECIFICATIONS

|                                                                   |                    |
|-------------------------------------------------------------------|--------------------|
| Model                                                             | BDIB-202           |
| Temperature control range                                         | RT. +5 °C~105 °C   |
| Timer range                                                       | 1 min-99 h 59 min  |
| Temp. stability @100 °C                                           | ≤0.5 °C            |
| Uniformity: within the block @37 °C / across similar block @37 °C | ±0.2 °C / ±0.3 °C  |
| Heating time (25 °C to 105 °C)                                    | ≤20 min            |
| Display accuracy                                                  | 0.1 °C             |
| Block quantity                                                    | 2                  |
| Dimension (LxWxH) mm                                              | 365x210x150 mm     |
| Weight                                                            | 4.5 kg             |
| Note                                                              |                    |
| Alt Name                                                          | Dry Bath Incubator |

## ACCESSORIES FOR PURCHASE

| No | Name  | Description                       |
|----|-------|-----------------------------------|
| 1  | Block | 96x0.2 mL PCR tube                |
| 2  | Block | 45x0.5 mL tube                    |
| 3  | Block | 35x1.5 mL tube                    |
| 4  | Block | 35x2.0 mL tube                    |
| 5  | Block | 15x0.5 mL+20x1.5 mL tube          |
| 6  | Block | 20x1.5 mL+15x2.0 mL tube          |
| 7  | Block | 32x0.2 mL+22x0.5 mL+9x1.5 mL tube |
| 8  | Block | 20x5 mL tube                      |
| 9  | Block | 20x10 mL tube                     |
| 10 | Block | 12x15 mL tube                     |
| 11 | Block | 6x50 mL tube                      |
| 12 | Block | 96x0.2 mL Elisa plate or 4 slides |
| 13 | Block | Solid block (no holes)            |
| 14 | Block | 384-well PCR plate                |
| 15 | Block | 40xΦ6 mm                          |

|    |       |                       |
|----|-------|-----------------------|
| 16 | Block | 28xΦ10 mm             |
| 17 | Block | 24xΦ12 mm             |
| 18 | Block | 24xΦ13 mm             |
| 19 | Block | 14xΦ15 mm             |
| 20 | Block | 14xΦ16 mm             |
| 21 | Block | 12xΦ19 mm             |
| 22 | Block | 11xΦ20 mm             |
| 23 | Block | 6xΦ26 mm              |
| 24 | Block | 6xΦ28 mm              |
| 25 | Block | 2xΦ40 mm              |
| 26 | Block | 2x7-12.5x12.5 cuvette |
| 27 | Block | Stainless plate       |



1



2



3



4



5



6



7



8



9



10




11



12



13



14



15



16



17



18



19



20



21



22



23



24



25



26



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)